

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SGK1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag2411

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MTVKTEAAKGLTYSRMRGMVAIIAFMKQRRMGLN  
DFIQKIANNYSYACKHPEVQSILKISQPQEPPELMNANP  
SPPSPSQQINLGPSSNPHAKPSDFHFLKVIKGSFGK  
VLLARHKAEVVFYAVKVLQKKAILKKKEEKHIMSERNV  
LLKNVKHPFLVGLHFSFQTADKLYFVLDYINGGELFYH  
LQRERCFLEPRARFYAAEIASALGYLHSLNIVYRDLKPE  
NILLDSQGHIVLTDGFLCKENIEHNSTSTFCGTPEYLA  
PEVLHKQPYDRTVDWWCLGAVLYEMLYGLPPFYSRN  
TAEMYDNI LNKPLQLKPNITNSARHLLGQLLQKDRKRL  
GAKDDFMEIKSHVFFSLINWDDLINKKITPPFNPNVSG  
PNDLRHFDPEFTEEPVNSIGKSPDSVLVTASVKEAAE  
AFLGFSYAPPTDSFL

(1-431aa encoded by BC001263)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.