

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human CYP4B1 protein



Basic Information

Catalog Number:

Ag2416

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SYYLAVSDLTLLMQQLVVSFYHNDFIYWLTPHGRRF
LRACQVAHDHTDQVIRERKAAALQDEKVRKKIQNRRH
LDFLDILLGARDEDDIKLSDADLRAEVDTFMFEGHDTT
TSGISWFLYCMALYPEHQHRCREEVREILGDQDFFQ
WDDLKGMTYLTMCIKESFRLYPPVPQVYRQLSKPVTF
VDGRSLPAGSLISMHIYALHRNSAVWPDPEVFDLRF
STENASKRHPFAFMPFSAGPRNCIGQQFAMSEMKVV
TAMCLLRFEFSLDPSRLPIKMPQLVLRSKNGFHLHLKPL
GPGSGK

(211-512 aa encoded by BC017758)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.