

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

MTHFD1 Fusion Protein



Basic Information

Catalog Number:

Ag24553

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

TNDRFLRKITIGQAPTEKGHTRTAQFDISVASEIMAVL
ALTTSLEDMRERLGKMMVVASSKKGEPVSAEDLGVSGA
LTVLMKDAIKPNLMQTLGTPVHVHAGPFANIAHGNS
SIIADQIALKLVGPEGFVVTEAGFGADIGMEKFFNIKC
RYSGLCPHVVVLVATVRALKMHGGGPTVTAGLPLPKA
YIQENLELVEKGFSNLKKQIENARMFGIPVVAVNAF
KTDTESELDLISRLSREHGAFDAVKCTHWAEGGKAL
ALAQAVQRAAQAPSSFQLLYDLKLPVEDKIRIIAQKIY
GADDIELLPEAQHKAEVYTKQGFNLPICTMAKTHLSL
SHNPEQKGVPTGFILPDIRASVGAGFLYPLVGTMTST
MPGLPTRPCFYDIDLDPETEQVNGLF
(536-935 aa encoded by BC009806)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.