

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SELL protein



## Basic Information

Catalog Number:

Ag2469

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGCRRTREGPSKAMIFPWKCQSTQRDLWNIFKLWG  
WTMLCCDFLAHHGTDCWTYHYSEKPMNWQRARRF  
CRDNYTDLVAIQNKAEIEYLEKTLPF SRSYYWIGIRKIG  
GIWTWVGTNKSLEEAENWGDGEPNNKKNKEDCVEI  
YIKRNKDA GKWND DACHKLKAALCYTASCQPWSCS  
GHGECVEIINNYTCNCDVGYGPPQCQFVIQCEPLEA  
PELGTMDCTHPLGNFSFSSQCAFSCSEGTNLGIEET  
TCGPFGNWSSPEPTCQVIQCEPLSAPDLGIMNCSHPL  
ASF SFTACTFICSEGTELGKKK TICESSGIWSNPSP  
CQKLDKSF SMIKEGDYNPLFIPVAVMV  
(1-355aa encoded by BC020758)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.