

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PIP₄K2A protein



Basic Information

Catalog Number:

Ag25139

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VMLMPDDFKAYSKI KVDNHLFNKENMPSHFKFKEYCP
MVFRLRERFGIDDQDFQNSLTRSAPLPNDSQARSQ
ARFHTSYDKRYIIKTITSEDVAEMHNILKKYHQYIVECH
GITLLPQFLGMYRLNVDGVEIYVIVTRNVFSHRLSVYR
KYDLKGSTVAREASDKEKAKELPTLKDNDFINEGQKIYI
DDNNKVFLEKLIKDVFLAQLKMDYSLLVGIHDVER
AEQEEVECEENDGEEEGESDGTHPVGTTPDSPGNTL
NSSPPLAPGEFDPNIDVYGIKCHENS PRKEVYFMAIIDI
LTHYDAKKAAHAAKTVKHGAGAEISTVNPEQYSKRF
LDFIGHILT

(59-406 aa encoded by BC018034)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.