

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human DUOXA1 protein



Basic Information

Catalog Number:

Ag2525

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MATLGHTFPFYAGPKTFPMDTTLASII MIFLTALATFI
VILPGIRGKTRLFWLLRVVTSLFIGAAILGTPVQQLNET
INYNEEFTWRLGENYAE EYAKALEKGLPDPVLYLAEKF
TPRSPCGLYRQYRLAGHYTSAMLWVAFLCWLLANVM
LSMPVLVYGGYMLLATGIFQLLALLFFSMATSLTSPCPL
HLGASVLHTHHGPAFWITLTTGLLCVLLGLAMAVHR
MQPHRLKAFFNQSVDEDPMLEWSPEEGLLSPRYRS
MADSPKSQDIPLSEASSTKAYCKEAHPKDPDCAL
(1-298aa encoded by BC020841)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.