

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# GABBR2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag25759

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET30a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MGLMPLTKEVAKGSI GRGVLPAVELAIEQIRNESLLRPY  
FLDLRLYDTECDNAKGLKAFYDAIKYGNHLMVFGGV  
CPSVTSIIAESLQGWNLVQLSFAATTPVLADKKKYPYF  
FRTVPSDNAVNPAILKLLKHYQWKRVGTLTQDVQRF5  
EVRNDLTGVLYGEDI EISDTEFSNDPCTSVKLLKGNL  
VRILLGQFDQNMAAKVFCAYEENMYGSIKQWIIPG  
WYEPSWWEQVHTEANSSRCLRKNLLAAMEGYIGVDF  
EPLSSKQIKTISGKTPQQYEREYNNKRSGVGPSKFHGY  
AYDGIWVIAKTLQRAMETLHASSRH  
(1 - 324 aa encoded by BC035071)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.