

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

STT3A Fusion Protein



Basic Information

Catalog Number:

Ag2698

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

STWVTSEAYSSPSIVLSARGGDSRIIFDDFREAYYWL
RHNTPEDAKVMSWWDYGYQITAMANRTILVDNNTW
NNTHISRVGQAMASTEELKAYEIMRELDVSYVLVIFGG
LTGYSSDDINKFLWMVRIGGSDTGKHIKENDYYTPTG
EFRVDREGSPVLLNCLMYKMCYYRFGQVYTEAKRPPG
FDRVRNAEIGNKDFELDVLEEAYTTEHWLVRIYKVKDL
DNRGLSRT

(475-705 aa encoded by BC020965)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.