

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PARG protein



Basic Information

Catalog Number:

Ag27288

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SNANIDRSPQNDDHSDTDSEENRDNQQLTTVKLAN
AKQTTEDEQAREAKSHQKCSKSCDPGEDCASCQQD
EIDVVPESPLSDVGSEVDGTGPKNDNLTRQESCLGN
SPPFEKESEPEPMDVDNSKNSCQDSEADEETSPGFD
EQEDGSSSQTANKPSRFQARDADIEFRKRYSTKGGEV
RLHFQFEGGESRTGMNDLNAKLPGNISSLNVECRNSK
QHGGKDSKITDHFMRPKAEDRRKEQWETKHQRTER
KIPK

(183-441 aa encoded by BC052966)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.