

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ARSD protein



## Basic Information

Catalog Number:

Ag27679

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ISWYSSFQFVRRWNCILMRNHDVTEQPMVLEKTASL  
MLKEAVSYIERHKHGPFLLFLSLHVIPLVTTSAFLGK  
SQHGLYGDNVEEMDWLIGKVLNAIEDNGLKNSTFTYF  
TSDHGGHLEARDGHSQ LGGWNGIYKGGKMGGWEG  
GIRVPGIFHWPGVLPAGRVIGEPTSLMDVFPVVQLV  
GGEVPPQDRVIDGHSLVPLLQGAEARSAHEFLFHYCG  
QHLHAARWHQKDSGSVWKVHYTTPQFHPEGAGAC  
YGRGVCPCSGEGVTHHRPPLFDLSRDPSEARPLTPDS  
EPLYHAVIARVGAAVSEHRQTLSPVPQQFMSNILW  
KPWLQPCCGHFPFCSCHEGDGTP  
(242-593 aa encoded by BC020229)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.