

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant mouse FC tag protein



## Basic Information

Catalog Number:

Ag27946

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PRGPTIKPCPPCKCPAPNLLGGPSVFI FPPKIKDVLMI S  
LSPIVTCVVVDVSEDDPDVQISW FVNNVEVHTAQT  
QTHREDYNSTLRVVSALPIQH QDWMSGKEFKCKVNN  
KDLPAPIERTISKPKGSVRA PQVYVLPPEEEMTKKQV  
TLTCMVTDFMPEDIYVEWTNNGKTELNYK NTEPVLDS  
DGSYFMYSKLRVEKKNWVERNSYSCSVVHEGLHNHH  
TTKFSRTPGK  
( encoded by )

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.