FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human DDR2 protein



Basic Information

Catalog Number:

Ag28183

Size:

50 μg

Form: Available lyophilized

Species: human

Expression Source:

 $\it e\,coli.$ -derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MFAKGVKIFKEVQCYFRSEASEWEPNAISFPLVLDDV NPSARFVTVPLHHRMASAIKCQYHFADTWMMFSEIT FQSDAAMYNNSEALPTSPMAPTTYDPMLKVDDSNT (291-398 aa encoded by NM_001014796)

Selected Validation Data