

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human IPAF protein



Basic Information

Catalog Number:

Ag30416

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EFRTLEVTLRDFS KLNKQDIRYLGKIFSSATSLRLQIKRC
AGVAGSLSLVLSTCKNIYSLMVEASPLTIEDERHITSVT
NLKTLSIHDLQNQR L PGGLTDSLGNLTKLIMDNIKM
NEEDA I KLA EGLK N L K M C L F H L T H L S D I G E G M D Y I V K
SLSSEPCDLEE I QLVSCCLSANAVKILAQNLHNLV KLSI
LDLSENYLEKDGNEALHELIDRMNVLEQLTALMLPWG
CDVQGSLSLLKHLEEV PQLVKLGLKNWRLTDEIRILG
AFFGKNPLKNFQQLNLAGNRVSSDGWLA FMGVFENL
KQLVFFDFSTKEFLPDPALVRKLSQVLSKLTFLQEARL
VGWQFDDDDLSVITGAFKLVTA

(658-1024 aa encoded by BC031555)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.