

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

spike protein Fusion Protein



Basic Information

Catalog Number:

Ag31496

Size:

50 µg

Form:

Available lyophilized

Species:

virus

Expression Source:

mammalian cells-derived, pHZ-k1sec, with N-terminal mychis.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RVQPTESIVRFPNITNLCPFGEVFNATRFASVYAWN
KRISNCVADYSVLYNSASFSTFKCYGVSPTKLNLCFT
NVYADSFVIRGDEVQRQIAPGQTGKIADYNYKLPDDFT
GCVIAWNSNNLDSKVGNYNYLRLFRKSNLKPFERDI
STEIYQAGSTPCNGVEGFNCYFPLQSYGFQPTNGVG
YQPYRVVVLSFELLHAPATVCGPKKSTNLVKNKCVNF
(319-541 aa encoded by NC_045512)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.