

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human 705CN protein



Basic Information

Catalog Number:

Ag31612

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MVAARALSVLLLLGLGGLVVADICPDDVDDPTDLCDO
EEDWSSGQRMAPNLNLLGLLAFVLPHSCLKHAALTA
AALCSLPGVAAVTCGDLKTFYKQQCCGSATKDLSA
TVPGCPYNFQPPACNMGEPTQPRDLSTNAEGDMTPK
AATLTDEQANFLPLVNVHFHLGAEHKSDSYNDGTD
EAYDAQSSGRRLAANPRPGWMCPTDTLTSSQLTSYN
FQYCEGVEVGKSYEVHYVHSSAGVDVDMTDGVNAD
DLADGLGGAANGRGILNPMVVVQGGVYQIVSGGPV
VDDLHGWTVVGHNSVMYPGSTTGQSWDNSVCS
PYAITWHVDKDC HQVSPESFDNLCKMMKEQYNMVS
DLKPHGSRILVDRKYVVKA EYVKPLA
(1-383 aa encoded by)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.