

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant alpaca VHH Mouse IgG2a Fc protein



Basic Information

Catalog Number:

Ag31696

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

1 QVQLVECGGG LVQPGGSLRL 21 SCVASGFTFS S
AYMSWVRQA 41 PGKGPEWVST ISTGGGIVNY 61 A
DSVKGRFAI SRDNAKNTLY 81 LQMNKLPED TALYY
CASNK 101 GPHYHSDYFD SNLYDFWGQG 121 TRV
TVSSGGC *

(1-130 aa encoded by)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.