

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PRKAG2 Fusion Protein



Basic Information

Catalog Number:

Ag3214

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EAGMLEKLEFEDEAVEDESEGVYMRFRMRSHKCYDIVP
TSSKLVVFDTTLQVKKAFFALVANGVRAAPLWESKK
QSFVGMILTDFINILHRYKSPMVQIYELEEHKIETWR
ELYLQETFKPLVNI SPDASLFDAVYSLIKNKIHRPVIDP
ISGNALYLTHKRILKFLQLFMSDMPKPAFMKQNLDEL
GIGTYHNI AF IHPDTP I I KALNIFVERRISALPVVDESGK
VVDIYSKFDVINLAAEKT YNNLDITVTQALQHR SQYF
EGVVKCNKLEI LETIVDRIVRAEVHRLVVVNEADSIVG
IISLSDILQALILTPAGAKQKETETE
(11-341 aa encoded by BC020540)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.