

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human AKT3 protein



## Basic Information

**Catalog Number:**

Ag32838

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MSDVTIVKEGWVQKRGEYIKNWRPRYFLKTDGSFIG  
YKEKPDVDLPYPLNFSVAKCQLMKTERPKPNTFIIR  
CLQWTTVIERTFHVDTPEREWEIAIQAVADRLQR  
QEEERMNCSPTSQIDNIGEEEMDASTTHHKRKT MND  
FDYLLKLGKTFGKVLVREKASGKYAMKILKKEVIAK  
DEVAHTLTERSVLKNTRHPFLTSKYSFQTKDRLCFVM  
EYVNGGELFFHLSRERVSEDRTRFYGAIVSALDYH  
SGKIVYRDLKLENLMRDKDGHIKITDFGLCKEGITDAAT  
MKTFCGTPEYLAPEVLEDNDYGRAVDWWGLGVV MY  
EMMCGRLPFYNQDHEKLFELILMEDIKFPRTLSSDAKS  
LLSGLLIKDPNKRLLGGPDDAKEIMRHSFFSGVNWQD  
VYDKLVPFPKQVTSSETDRYFDEEFTAQTITITPPE  
KCQQSDCGMLGNWKK  
(1-465 aa encoded by NM\_005465)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.