

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

STK33 Fusion Protein



Basic Information

Catalog Number:

Ag3561

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MADSGLDKKSTKPCDCSSASQKDVLCVCSKTRVPP
VLVVEMSQTSSIGSAESLISLERKKEKNINRDITSRKDL
PSRTSNVERKASQQQWGRGNFTEGKVPHIRIENGAA
IEEIYTFGRILGKGSFGIVIEATDKETETKWAIKV NKEK
AGSSAVKLLEREVNILKSVKHEHIIHLEQVFETPKKMYL
VMELCEDGELKEILDRKGFSENETRWIIQSLASAIAYL
HNNDIVHRDLKENIMVKSSLI DDNNEINLNKVTDFGL
AVKKQSRSEAMLQATCGTPIYMAPEVISAHDYSQQC
DIWSIGVVMYMLLRGEPPLASSEEKLFELIRK GELHFE
NAVWNISIDCAKSVLKQLMKVDPAHRITAKELLDNQ
WL

(1-381 aa encoded by BC031231)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.