

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SIRT7 Fusion Protein



## Basic Information

**Catalog Number:**

Ag3656

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAAGLSRSEKAAERVRLREEQQRERLRQVSRLRK  
AAASAEAEGRLLAESADLVTELQGRSRRREGLKRRQE  
EVCDDPEELRGKVELASAVRNAKYLVVYTGAGISTA  
ASIPDYRGPNGVWTLQKGRSVSAADLSEAEPTLTHM  
SITRLHEQKLVQHVVSNQCDGLHLRSLPRTAISELHG  
NMYIEVCTSCVPNREYRVFQVTERALHRHQGTGR  
CHKCGTQLRDTIVHFGERTLGQPLNWEAAEAAASR  
ADTILCLGSSLVKLYPRLWCMKPPSRPKLYIVNLQ  
WTPKDDWAALKLHGKDDVMRLMAELGLEIPAYSR  
WQDPFSLATPLRAGEEGSHSRKSLCRSREEAPPDGR  
GAPLSSAPILGGWFRGCTKRTKRKKT  
(1-400 aa encoded by BC017305)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.