

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PIAS4 protein



Basic Information

Catalog Number:

Ag37383

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PNIAVKVNHSYCSVPGYYPSNKPGEVKRPCRPINLTH
LMYLSSATNRITVTWGNYGKSYVALYLVRQLTSSELL
QRLKTIGVKHPELCKALVKEKLRDPDSEIATTGVRVSL
ICPLVKMRLSVPCRAETCAHLQCFDAVFYLMNEKKP
TWMCPVCDKPAPYDQLIIDGLLSKILSECEDADEIYL
VDGSWCPIRAEKERSCSPOGAILVLGPSDANGLLPAP
SVNGSGALGSTGGGPGVSMENKPGADVVDLTLDLS
SSSSEDEEEEEEEDEEGPRPKRRCPFQKGLVPAC
(210-510 aa encoded by BC029874)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.