

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ASNSD1 protein



Basic Information

Catalog Number:

Ag37503

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SVVANEAKLYLEKPVVPLNMMLPQAALETHCSNISNV
PPTREILQVFLTDVHMKEVIQQFIDVLSVAVKKRVLCL
PRDENLTANEVLKTCDRKANVAIFSGGIDSMVIATLA
DRHIPLDEPIDLLNVAFIAEETMPTTFNREGNKQKNK
CEIPSEEFKDVAAAAADSPNKHVSVPDITGRAGLK
ELQAVSPSRIWNFVEINVSMEEQLKRRRTRICHLIRPLD
TVLDDSIGCAVWFASRGIGWLVAQEGVKSYSNAK
VVLTGIGADEQLAGYSRHRVRFQSHGLEGLNKEIMME
LGRISRNLRDDDRVI GDHGKEARFPFLDENNVVSFLNS
LPIWEKANLTPRGI GEKLLRLAAVELGLTASALLPKRA
MQFGSRIAKMEKINEKASDKCGRQLIMSLNLSIEKET
KL

(227-643 aa encoded by BC001243)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.