

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human NTRK2 protein



Basic Information

Catalog Number:

Ag3790

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SRIWCDSPSPGIVAFPRLEPNSVDPENITEIFIANQKRL
EII NEDDVEAYVGLRNLTIVDSGLKFVAHKAFKNSNL
QHINFTRNKLTSLSRKHFRLDLSEILVGNPFTCSCDI
MWIKTLQEAKSSPDTQDLYCLNESSKNIPLANLQIPNC
GLPSANLAAPNLVVEEGKSITLSCSVA G D P V P N M Y W D
VGNLVS K H M N E T S H T Q G S L R I T N I S S D D S G K Q I S C V A
ENLVGEDQDSVNLTVHFAPTITFLESPTS D H H W C I P F
TVKGNPKPALQWFYNGAILNESKYICTKIHVTNHTEY
HGCLQLDNPTHMNGDYTLIAKNEYGKDEKQISAHF
MGWPGIDDGANPNYPDVIYEDYGTAAANDIGDTTNR
NEIPSTDVTDKGTGREHLSV

(41-433 aa encoded by BC031835)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.