

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PDS5B Fusion Protein



## Basic Information

Catalog Number:

Ag3800

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

AI EPYITNFFNQV LMLGKTSISDLSEHVFDLILELYNIDS  
HLLSVLPQLEFKLSNDNEERLQVVKLLAKMFGAKDS  
ELASQNKPLWQCYLGRFNDINVPIRLEC VKFASHCLM  
NHPDLAKDLTEYLVKVRSHDPEEAI RHDVIVSIVTAAKK  
DILLVNDHLLNFVRERTLDKRWRVRKEAMMGLAQIYK  
KYALQSAAGKDAAKQIAWIKDKLLHIYYQNSIDDRLLV  
ERIFAQYMPVPHNLETTERMKCLYYLYATLDLNAVKYVS  
NIKFCSFHPLQYIGFYGKETTNTCILKCNLCNVNIV  
(228-529 aa encoded by BC039256)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.