

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human TAF1B protein



## Basic Information

Catalog Number:

Ag38211

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MDLEESEEFKERTQCAAVSWGLTDEGKYYCTSCHN  
VTERYQEVNTNDLIPNTQIKALNRGLKKKNTEKGWD  
WYVCEGFQYILYQQAELKNLGCVGPELKNLHNFV  
KRYLQKSKQAYCKNPVYTTGRKPTVLEDNLSHSDWAS  
EPELSDVSCPPFLESQSDIHTRKPFVSKASQS  
ETSVCSGLDGVSEYSQRKEKGIKMTMPQTLAFCYLS  
LLWQREAITLSDLLRFVEEDHIPYINAFQHFPEQMKLY  
GRDRGIFGIESWPDYEDIYKKTIEVGTFLDLPRFPDITE  
DCYLHPNILCMKYLMEVNLPEMHSLTCHVVKMTGM  
GEVDFLTFDPIA

(1-346 aa encoded by BC018137)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.