

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RAD51L1 Fusion Protein



Basic Information

Catalog Number:

Ag3884

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MGSKLKRVLGSLQELCDRLSRHQILTCQDFLCLSPLEL
MKVTGLSVYRGVHELLCMVSRACAPKMQTAYGIKAQR
SADFSPAFLSTTLSALDEALHGGVACGSLTEITGPPGC
GKTQFCIMMSILATLPTNMGLEGAVVYIDTESAFSA
ERLVEIAESRFPRYFNTEEKLLTSSKVHLYRELTCDLVL
QRILEEEIISKGIKLVILDSVASVVRKEFDAQLQGNL
KERNKFLAREASSLKYLAEFESIPVILTNIITTHLSGALA
SQADLVSPADDLSLSEGTSGSSCVIAALGNTWSHSV
NTRLILQYLDERRQILIAKSPAPFTSFVYTIKEEGLVL
QETTFCSVTQAEINWAPEILPPQPPEQLGLQMCHHT
QLIF

(1-384aa encoded by BC030219)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.