

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Galc Fusion Protein



Basic Information

Catalog Number:

Ag3914

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

YYLKTVGHLEKGGSYVALTDGLGNLTIIIETMSHQHSM
CIRPYLPYYNVSHQLATFTLKGSLREIQELQVWYTKLG
TPQQRHFQQLDTLWLLDGGSGFTLELEDEIFTLTLT
TGRKGSYPPPPSSKPFPTNYKDDFNVEYPLFSEAPNFA
DQTGVFEYMNEDREHRFTL

(340-513 aa encoded by BC086671)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.