

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FOXI1 protein



Basic Information

Catalog Number:

Ag3958

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSSFDLPA^{SP}PRCSPQ^FPSIGQEPPEMNL^YYENFFH
PQGV^{PS}QRP^SFE^{GG}GEYGATPN^PYLW^FNGPTMT^PP
PYLP^{GP}NAS^PFLPQAYGV^QRPL^{PS}VSG^LGGSD^LGW^L
PIPSQEE^LMK^LVRPP^SYSALIA^{MA}IHGAP^DKRL^TLSQ^I
YQYV^{AD}NFP^FY^NKS^KAGWQ^NSIRH^NLSL^NDC^FKK^VP
RDE^{DD}PA^YVSG^SPTS^HPL^VTP^GLSPE^SDK^TGQ^NSL
TFNS^FSPL^TNLS^NHSG^GGDW^ANPM^PTN^MLSY^GGS^VL
SQ^FSP^HFY^NSV^NTS^GVLY^PREG^TEV
(1-283 aa encoded by BC029778)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.