

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MIA3 protein



## Basic Information

Catalog Number:

Ag39740

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KAVSAAEEVKTYKRRIEMEDELQKTERSFKNQIATHE  
KKAHENWLKARAAERAIAEEKREANLRHKLELTQK  
MAMLQEEPVIKMPGKPNQNPVRRCPVLSQNGSFG  
PSPVSGGECSPPLTVEPPVRPLSATLNRRDMPRSEFGS  
VDGPHPRWSAEASGEPSPDPGSGTAMMNSSSR  
GSSPTRVLEDEKVNMAPKGGPPFPGVPLMSTPMGGP  
VPPPIRYGPPQVLCGPFGRPLPPFPGMRPPLGLRE  
FAPGVPPGRRDPLHPRGFLPGHAPFRPLGSLGPREYFI  
PGTRLPPTHGQVEYPPPAVRDLLPSGSRDEPPASQ  
STSQDCSQALKQSP

(118-467 aa encoded by BC047116)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.