

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

BPIL1 Fusion Protein



Basic Information

Catalog Number:

Ag4051

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

II LPTD A T P F V L P R H V G T E G S M A T G G L S Q Q L F D S A L L L
L Q K A G A L N L D I T G Q L R S D D N L L N T S A L G R L I P E V A R Q F
P E P M P V V L K V R L G A T P V A M L H T N N A T L R L Q P F V E V L A
T A S N S A F Q S L F S L D V V V N L R L Q L S V S K V K L Q G T T S V L
G D V Q L T V A S S N V G F I D T D Q V R T L M G T V F E K P L L D H L N
A L L A M

(233-424 aa encoded by BC034415)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.