

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TDG Fusion Protein



Basic Information

Catalog Number:

Ag4190

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEAENAGSYSLQQAQAFYTFPFQQLMAEAPNMAV
VNEQQMPPEVVPAPAPAQEPVQEAPKGRKRKPRTEP
KQPVEPKKPVESKKGSAKSEKQEKITDTFKVKRKV
DRFNGVSEAEELTKLPDILTFNLDIVIIGINPGLMAAYK
GHHYPGPGNHFWKCLFMSGLSEVQLNHMDDHTLPG
KYGIGFTNMVERTTPGSKDLSSKEFREGGRILVQKLQK
YQPRIAVFNGKCIYEIFSKEVFGVKVKNLEFGLQPHKI
PDTE TLCYGPSSSARCAQFPRAQDKVHYYIKLKDRL
DQLKGIERNMDVQEVQYTFDLQLAQEDAKKMAVKE
EKYDPGYEAA YGGAYGENPCSSEPCGFSSNGLIESVE
LRGESAFSGIPNGQWMTQSFTDQIPSFNHC GTQE Q
EEESHA

(1-410 aa encoded by BC037557)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.