

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human TDRD3 protein



## Basic Information

**Catalog Number:**  
Ag4351

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
YRSSNTEQNGVKDNNHLRHPPRNDTRQPRNEKPPRF  
QRDSQNSKSVLEGSGLPRNRGSRPSTSSVSEVWAE  
DRIKCDRPSRYDRTKDTSYPLGSQHSDFKRDNS  
MQSRSGKGPSFAEAKENPLPQGSVDYNNQKRGKRES  
QTSIPDYFYDRKSQTINNEAFSGIKIEKHFVNTDYQN  
PVRNSNFIGVNPGEVEMPLKGRRI GPIKAGPVTAVPC  
DDKIFYN SGPKRRSGPIKPEKILESSIPMEYAKMWKPG  
DEC FALY WEDNKFYRAEVEALHSSGMTAVVKFIDYG  
NYEEVLLSNIKPIQTEAWEEEGTYDQTLFRRGGDGO  
PRRSTRPTQQFYQPPRARN  
(301-651 aa encoded by BC030514)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.