

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human THADA protein



Basic Information

Catalog Number:

Ag4414

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

TSSALSMGPFVPFIMRCGHSPVYHSREMAARALVPFV
MIDHIPNTIRSLSTLPSC TDQCFRQNHIGTLLQVFH
LLQAYSDSKHGTNSDFQHELT DITVCTKAKLWLA KRQ
NPCLVTRAVYIDIFLLTCLNRS AKDNQPVLES LGFW
EEVRGII SGSELTGFPWAFKVPGLPQYLSLRLAIAA
VWAAAAKSGERETNPISFSQLLES AFPEVRS LLEAL
LEKF LAAA SGLGEKGV PLLCNMGEK FLLAMKENHPE
CFCKILKILH C MDPGEWLPQTEHC VHLTPKEFLI WTM
DIASNERSEIQSVALRLASKVISHHMQTCVENRELI AA
EL

(59-400aa encoded by BC025773)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.