

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SENP2 Fusion Protein



Basic Information

Catalog Number:
Ag4664

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MKETIQGTGSWGPEPPGPIPPAYSSPRRERLRWPPP
PKPRLKSGGGFGPDGSGTTVPARRLPVPRPSFDASA
SEEEEEEEEEDEEEEEVAAWRLPPRWSQLGTSQRP
RPSRPTHRTCSQRRRRAMRAFRMLLYSKSTSLTFHW
KLWGRHRGRRRGLAHPKNHLSQGGATPQVPSGCC
RFDSPRGPPPRLLGALMAEDGVRGSPVPVPSGPPM
EEDGLRWTPKSPDPDSELLSCTLPNGFGGQSGPEGE
RSLAPPDASILISNVCSIGDHVAQELFQGSDLGMAEE
AERPGEKAGQHSPLREEHVTCVQSILDEFLOQTYGSLIP
LSTDEVVEKLEDIFQQEFSTPSRKLGLVLIQSYQRMP
GNAMVRFGRVAYKRHVLTMDLGLTYGQNLWLDQV
MNMVYDGLVMDTVPEKVHFFNSFFYDKLRTKGYDGVK
RWTKNVDIFNKELLLPIHLEVHWSLISVDVRRRTITYF
DSQRTLNRRCPIKHIKYLQAEAVKDKRDLDFHQGWKG
YFKMNVARQNNDSDCGAFVLQYCKHLALSQPF5FT
QQDMPKLRQIYKELCHCKLTV
(1-589aa encoded by BC080658)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.