

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SETX Fusion Protein



Basic Information

Catalog Number:

Ag4785

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LVMEIHKLIKDKRKDVSRNIGIITHYKAQKTMIQKDL
KEFDRKGPAEVDTVDAFQGRQKDCVIVTCVRANSIQ
GSI GFLASLQRLNVTITRAKYSLFILGHLRTL MENQHW
NQLIQDAQKRGAIKTCDKNYRHDAVKILKLPVLQRS
LTHPPTIAPEGSRPQGGLPSSKLD SGFAKTSVAASLYH
TPSDSKEITLTVTSKDPERPPVHDQLQDPRLLKRMGIE
VKGGIFLWDPQPSSPQHPGATPPTGEPGFVHVHQL
SHIQQPAAVVAALSSHKPPVRGEPAA SPEASTCQS
KCDDPEEELCHRREARAFSEGEQEKCGSETHHTRRNS
RWDKRTLEQEDSSSKRKLL

(575-930 aa encoded by BC032600)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.