

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# DARS2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag4809

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

LI E G L L Q Y S W P N D K D P V V V P F P T M T F A E V L A T Y G T D K  
P D T R F G M K I I D I S D V F R N T E I G F L Q D A L S K P H G T V K A I  
C I P E G A K Y L K R K D I E S I R N F A A D H F N Q E I L P V F L N A N R  
N W N S P V A N F I M E S Q R L E L I R L M E T Q E E D V V L L T A G E H  
N K A C S L L G K L R L E C A D L L E T R G V V L R D P T L F S F L W V V D  
F P L F L P K E E N P R E L E S A H H P F T A P H P S D I H L L Y T E P K K A  
R S Q H Y D L V L N G N E I G G G S I R I H N A E L Q R Y I L A T L L K E D V  
K M L S H L L Q A L D Y G A P P H G G I A L G L D R L I C L V T G S P S I R D  
V I A F P K S F R G H D L M S N T P D S V P P E E L K P Y H I R V S K P T D  
S K A E R A H

(296-645 aa encoded by BC045173)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.