

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

MAP2K3 Fusion Protein



Basic Information

Catalog Number:

Ag4861

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MESPASSQPASMPQSKGKSKRKKDLRISCMKPPAPN
PTPPRNLDsrTFITIGDRNFVEADDLVTISELGRGAYG
VVEKVRHAQSGTIMAVKRIRATVNSQEQKRLMLDLI
NMRTVDCFYTVTFYGALFREGDVVICMELMDTSLDK
FYRKVLDKNMTIPEDILGEIAVSIVRALEHLHSLSVIHR
DVKPSNVLIINKEGHVKMCDFGISGYLVDSVAKTMDA
GCKPYMAPERINPELNQKGYNVKSDVWSLGITMIEMA
ILRFPYESWGT P F Q L K Q V V E E P S Q L P A D R F S P E F V
D F T A Q C L R K N P A E R M S Y L E M E H P F F L H K T K I K T D I A A
FVKEILGEDS

(1-347 aa encoded by BC032478)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.