

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BTBD1 protein



Basic Information

Catalog Number:
Ag4893

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MASLGPAAAGEQASGAEAPGAPPPPPSPSLGPL
LPLQREPLYNWQATKASLKERFAFLNSELLSDVRFVL
GKGRGAAAAGGPQRIPAHRFVLAAGSAVFDAMFNG
GMATTSAEIELPDVEPAFLALLRFLYSDEVQIGPETV
MTTLYTAKKYAVPALEAHCVEFLT KHLRADNAFMLLT
QARLFDEPQLASLCLDTIDKSTMDAISAEGFTDIDIT
LCAVLERDLSIRESRLFVAVRWAEECQRQLPVT
FGNKQKVLGKALSIRFPLMTIEEFAAGPAQSGILSDR
EVVNLFHFTVNPKPRVEYIDRPRCCLRGKECCINRFQ
QVESRWGYSGSDRIRFTVNRRI SIVGFGLYGSIHGPT
DYQVNIQIIEYEKKQLTGQNDTGFSCDGTANTFRVM
FKEPIELPNVCYTACATLKGPD SHYGTGKGLKVVHET
PAASKTVFFFFSSPGNNNGTSIEDGQIPEIIFYT
(1-484aa encoded by BC028097)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.