

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SLC12A2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag4991

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

GSSTQALTYLNALQHSIRLSGVEDHVKNFRPQCLVMT  
GAPNSRPALLHLVHDFTKNVGLMICGHVHMGP RRQA  
MKEMSIDQAKYQRWLIKKNMKAFYAPVHADDLREGA  
QYLMQAAAGLGRMKPNTLVLGFKDWLQADM RDVDM  
YINLFHDAFDIQYGVVIRLKEGLDISHLQGG EELLSS  
QEKSPGTDKDVVVSVEYSKSDLDTSKPLSEK PITHKVE  
EEDGKTATQPLLKESKGPVPLNVADQKLEASTQF  
QKKQGKNTIDVWWLFDGGLLIPYLLTTKKKWKDC  
KIRVFIGGKINRIDHRRAMATLLSKFRIDFSDIMVLGD  
INTKPKKENIAFEEIIEPYRLHEDDKEQDIADKMK EDE  
PWRTDNELELYKTKFYEP C  
(759-1150 aa encoded by BC033003)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.