

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# LRRC4C Fusion Protein



## Basic Information

**Catalog Number:**

Ag5035

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

QTCPSVCSCSNQFSKVICVRKNLREVPDGISTNTRLL  
NLHENQIQIHKVNSFKHLRHLEILQLSRNHIRTIEIGAFN  
GLANLNTLELFDNRLTTPNGAFVYLSKLELWLRNNPI  
ESIPSYAFNRIPSLRRLDLGELKRLSYISEGAFEGLSNLR  
LNLAMCNLREIPNLTPLIKLELDLSGNHLSAIRPGSFQ  
GLMHLQKLWMIQSQIQVIERNAFDNLQSLVEINLAHN  
NLTLPHDLFTPLHHLERIHLLHNPWNCNCDILWLSW  
WIKDMAPSNTACCARCNTPPNLKGRYIGELDQNYFT  
C

(45-351 aa encoded by BC041374)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.