

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SLC22A6 protein



## Basic Information

Catalog Number:

Ag5072

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

TAAIPTHHCRPPADANLSKNGGLEWLPDRDQGQPE  
SCLRFTSPQWGLPFLNGTEANGTGATEPCTDGIWYD  
NSTFPSTIVTEWDLVCSHRALRQLAQSLYMGVLLGA  
MVFGYLADRLGRRKVLILNYLQTAVSGTCAAFAPNFPI  
YCAFRLLSGMALAGISLNCMTLNVEWMPHTRACVGT  
LIGVYVSLGQFLLAGVAYAVPHWRHLQLLSAPFFAF  
FIYSWFFIESARWHSSGRDLTLRALQRVARINGKRE  
EGAKLSMEVLRASLQKELTMGKGQASAMELLRCPTLR  
HLFLC

(40-340aa encoded by BC033682)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.