

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human GRM3 protein



## Basic Information

Catalog Number:

Ag5306

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HNFLRREIKIEGDLVLGGFLPINEKGTGTEECGRINEDR  
GIQRLEAMLF AIDEINKDDYLLPGVKLGVHILDTC SRDT  
YALEQSLEFVRA SLTKVDEAEYMC PDGSYAIQENIPLL  
IAGVIGGSYSSVSIQVANLLRLFQIPQISYASTSAKLSL  
KSRDYDFARTVPPDFYQAKAMAEILRFFNWTYVSTVA  
SEG DYGETGIEAFEQEARLRNICIATAEKVGRSNIRKS  
YDSVIRELLQKPNARVVVLFMRSDDSRELI AAA SRAN  
ASFTWVASD GWGAQESI IKGSEHVAYGAITLEASQ  
PVRQFD RYFQSLNPNYNNHRNPWFRDFWEQKFC SL  
QNKRNHRRVCDKHLAIDSSNYEQESKIMFVVNA  
(26-396 aa encoded by BC022496)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.