

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MORC1 protein



## Basic Information

Catalog Number:

Ag5396

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EETDSDVEYISETKIMKKSMECKMNSQQQRIPVALPE  
NVKLAERSQRSQIANITTVWRAQPTEGCLKNAQAAS  
WEMKRRKQSLNFVEECKVLTEDENTSDSDIILVSDKSN  
TDVSLKQEKKEIPLLQEKQELCNDVLA MKRSSLPS  
WKSLLNVP MEDVNLSSGHIA RVSVSGCKV ASSPASS  
QSTPVKETVRKLSKLR EILLYFFPEHQLPSELEEPALSC  
ELEQCPEQM NKKLK MCFNQIQNTY MVQYEKKIKRKL  
QSIIYDSNTRGIHNEISLGQCENKRKISEDKLNLR IKLA  
LLLQKLQLGGPEGDLQTDTYLEALLKEDNLLFQNNLN  
KVTIDARHRLPLEKNEKTS EN  
(626-984aa encoded by BC050307)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.