

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TBC1D1 Fusion Protein



Basic Information

Catalog Number:

Ag5406

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LAEQFHLKHQFPSKQQPKDVVPYKELLKQLTSQQHAIL
IDLGRTFPTHYPYFSAQLGAGQLSLYNIKAYSLLDQEV
GYCQGLSFVAGILLHMSEEEAFKMLKFLMFDMLRK
QYRPDMIIQIQMYQLSRLLDYHRDLYNHLEEHEIGP
SLYAAPWFLTMFASQFPLGFVARVFDMIFLQGTVEVIF
KVALSLLGSHKPLILQHENLETIVDFIKSTLPNLGLVQM
EKTINQVFEMDIKQLQAYEVEYHVLQEEELIDSSPLSD
NQRMDKLEKTNSSLRKQNLDLLEQLQVANGRIQSLEA
TIEKLLSSEKLLQAMLTLELERSALLQTV EELRRRSAEP
SDREPECTQPEPTGD

(813-1168 aa encoded by BC050321)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.