

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

WASL Fusion Protein



Basic Information

Catalog Number:

Ag5424

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSSVQQQPPPPRRVTNVGSLLLTPQENESLFTFLGKK
CVTMSSAVVQLYAADRNCMWSKKCSGVA CLVKDNP
QRSYFLRFIDIKDGKLLWEQELYNNFVYNSPRGYFHTF
AGDTCQVALNFANEEAAKFRKAVTDLLGRRQRKSEK
RRDPPNGPNLPMATVDIKNPEITTRFYGPQVNNISH
TKEKKKGKAKKRLTKADIGTPSNFQHIGHVGDWPNP
GFDLNNLDPELKNLDMCGISEAQLKDRETSKVIYDFI
EKTGGVEAVKNELRRA P P P P P S R G P P P P P P P H N
S G P P P P P A R G R G A P P P P P S R A P T A A P P P P P S R P S V A
V P P P P P N R M Y P P P P P A L

(1-350 aa encoded by BC052955)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.