

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

NEK2 Fusion Protein



Basic Information

Catalog Number:

Ag5474

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MPSRAEDYEVLTYTIGTGSYGRCQKIRRKSDGKILVWKE
LDYGSMTAEAEKQMLVSEVNLLRELKHPNIVRYDRIID
RTNTTLYIVMEYCEGGDLASVITKGTKERQYLDEEFVL
RVMTQLTLALKECHRRSDGGHTVLHRDLKPANVFLDG
KQNVKLGDFGLARILNHDTSF AKTFVGTPTYMSPEQM
NRMSYNEKSDIWSLGCLLYELCALMPPFTAFSQKELA
GKIREGKFRRI PYRYSDELNEIITRMLNLKDYHRPSVEEI
LENPLIADLVADEQRRNLERRGRQLGEPKSDSSPVL
SELKLEIQLQERERALKAREERLEQKEQELCVRERLAE
DKLARAENLLKNYSLLKERK

(1-362 aa encoded by BC043502)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.