

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ZHX1 protein



Basic Information

Catalog Number:

Ag5770

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SEQDPDELI S D L D E G P P V L T P V E N T R A E S I S S D E E V H
E S V D S D N Q Q N K K V E G G Y E C K Y C T F Q T P D L N M F T F H
V D S E H P N V V L N S S Y V C V E C N F L T K R Y D A L S E H N L K Y H
P G E E N F K L T M V K R N N Q T I F E Q T I N D L T F D G S F V K E E N
A E Q A E S T E V S S S G I S I S K T P I M K M M K N K V E N K R I A V H
H N S V E D V P E E K E N E I K P D R E E I V E N P S S S A S E S N T S T S I
V N R I H P S T A S T V V T P A A V L P G L A Q M I T A V S A Q Q N S N
L I P K V L I P V N S I P T Y N A A L D N N P L L N T Y N K F P
(16-307 aa encoded by BC040481)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.