

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TCF7 Fusion Protein



## Basic Information

**Catalog Number:**

Ag5989

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MPQLDSGGGGAGGGDDLGAPELLAFQDEGEEDD  
KSRDSAAGPERDLAELKSSLVNESEGAAGGAGIPGVP  
GAGAGARGEAEALGREHAAQRLFPDKLPEPLEDGLKA  
PECTSGMYKETVYSAFNLLMHYPPPSGAGQHPQPQP  
PLHKANQPPHGVPLSLYEHFNSPHPTAPADISQKQ  
VHRPLQTPDLSGFYSLTSGSMGQLPHTVSWFTHPSLM  
LGSVPGHPAAIPHPAIVPPSGKQELQPFDRNLKTA  
ESKAEKEAKKPTIKKPLNAFMLYMKEMRAKVIAECLK  
ESAAINQILGRRWHALSREEQAKYYELARKERQLHMQ  
LYPGWSARDNYGKKKRRSREKHQESTTGGKRNAFGTY  
PEKAAAPAPFLPMTVL

(1 - 384 aa encoded by BC048769)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.