

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# DYX1C1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag5997

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MPLQVSDYSWQQTKTAVFLSLPLKGVCRDTRDVFCT  
ENYLKVNFPPLFEAFLYAPIIDDESSKAKIGNDTIVFTL  
YKKEAAMWETLSVTGVDKEMMQRIREKSI LQAQERA  
KEATEAKAAAKREDQKYALSVMMKIEEEERKKIEDMK  
ENERIKATKALEAWKEYQRKAEQKKIQREEKLCQKE  
KQIKEGRKIKIKYSLTRNLASRNLA PKGRNSENIFTEKLIK  
EDSIPAPRSVGSIKINFTRVFPALRESQVAEEEEWL  
HKQAEARRAMNTDIAELCDLKEEKNPEWLKDKGNKL  
FATENYLAAINAYNLAIRLNNKMPLLYNRAVCHLKLK  
NLHKAIEDSSKEFC SLEGIECQASEPKLSHHIPSDLHVY  
IQMA

(1-381 aa encoded by BC062564)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.