

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TOP3A Fusion Protein



## Basic Information

**Catalog Number:**

Ag6010

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

PEPIRKCPQC NKDMVLKTKKNGGFYLS CMGFPECRSA  
VWLPDSVLEASRDSSVCPVCQPHPVYRLKLFKRGSL  
PPTMPLEFVCCI GGDDTLREILD LRFSGGPPRASQPS  
GRLQANQSLNRMDNSQH PQPADSRQTGSSKALAQ T  
LPPPTAAGESNSVTCNC GQEA VLLTVRKEGPNRGRQ  
FFKCNNGGSCNFFLWADSPNPGAGGPPALAYRPLGAS  
LGCPPGPIHLGGFGNPGDGS SGTSC LCSQPSVTRT  
VQKDGPNGRQFHTCAKPREQQCGFFQW DENTA  
PGTSGAPSWTGDRGR TLESEARSKRPRASSSDMGSTA  
KKPRKCSLCHQPGHTRPFCPQNR

(652-1001 aa encoded by BC051748)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.